MAGUS CHD10 DIGITAL CAMERA

MAGUS

www.magusmicro.com

Levenhuk Inc. (USA):
928 E 124th Ave. Ste D, Tampa, FL 33612, USA,
+1-813-468-3001,
contact_us@levenhuk.com

Levenhuk Optics s.r.o. (Europe):
V Chotejně 700/7, 102 00 Prague 102, Czech Republic,
+420 737-004-919,
sales-info@levenhuk.cz

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EN MAGUS CHD10 DIGITAL CAMERA

Never look at the sun, a bright light source, or laser beam through the camera. It can cause permanent damage to your eyes.

MAGUS high-speed HDMI cameras are designed to be paired with a digital microscope. Sony Starvis high-sensitivity sensor ensures perfect color reproduction and image quality even in low light conditions. You can install the camera software to view, save, and edit captured images and videos. The rugged aluminum housing is designed for operation at low temperatures (-10°C).

CAMERA INSTALLATION

1. Use an HDMI cable to connect the camera to an HDMI monitor (1).
2. Connect a USB mouse to the connector (2).
3. Insert the memory card into a memory card slot (3).
4. Connect the camera to the power supply (4).
5. Turn the monitor on and watch videos in the built-in program. Use your mouse to control the camera.

CALIBRATION USING A CALIBRATION SLIDE


You should use a special stage micrometer (calibration slide) to determine the linear dimensions (in millimeters or microns).

The calibration slide is a transparent glass (of the same size as the specimen slide) that has a micrometer scale with a scale division of 0.01mm etched on the surface.

The calibration slide (stage micrometer) is used to calibrate the image analysis software for measurements in actual units. In the calibration mode, you should capture an image of the micrometer scale with every objective magnification and indicate the known distance. That lets you establish a scale of the image in actual units (micrometer, millimeter, etc.).

Calibration:


1. Place the calibration slide on the microscope stage.
2. Select the required objective and set the maximum camera resolution.
3. Get a contrast image of the scale on the monitor screen and capture the image.
4. Select the "Calibrate" function in the software you are using.
5. Double-click on the maximum visible distance and enter the value in actual units.
6. Enter the calibration setting and check the result. The program will save the calibration factor.
7. You can select any measurement unit later, and all the results will be re-calculated in accordance with this selection.

CONTROL PANEL

To open the Camera control panel, move the mouse cursor to the left of the window. The control panel will automatically pop up.

CAMERA CONTROL PANEL

FUNCTION

FUNCTION DESCRIPTION

Snap	Captures the image and saves it to the SD card.
Record	Records a video and saves it to the SD card.
Auto Exposure	Automatically adjusts the exposure time and gain according to the exposure compensation value.
Exposure Compensation	Available when Auto Exposure is selected. Adjusts the exposure compensation according to the current video brightness to achieve the proper brightness value.
Exposure Time	Available when Auto Exposure is not selected. Adjusts the time of exposure and thus adjusts the brightness of the video.
Gain	Adjusts the brightness of the video. The Noise will be adjusted accordingly.
Red	Adjusts the proportion of the red color in RGB in the video.
Green	Adjusts the proportion of the green color in RGB in the video.
Blue	Adjusts the proportion of the blue color in RGB in the video.
Auto White Balance	Adjusts the White Balance of the video continuously.
Manual White Balance	Adjusts the Red color or Blue color parameters to set the White Balance of the video. White balance is adjusted to obtain a natural ratio of cold and warm colors on the image with regard to the light source type. Move the corresponding sliders to adjust the values.
ROI White Balance	Adjusts the white balance of the Region of Interest when the Region of Interest content is changed.
Sharpness	Adjusts the Sharpness level of the video.
Denoise	Removes the noise from the video.
Saturation	Adjusts the Saturation level of the video.
Gamma	Allows you to adjust the midtone values of the image. The higher the Gamma value is, the darker the image will be.
Contrast	Adjusts the Contrast level of the video.
DC (Direct current)	Allows you to use a light source without light flickering.
AC(50HZ) (Alternative current)	Eliminates light flickering when a 50Hz light source is used.
AC(60HZ) (Alternative current)	Eliminates light flickering when a 60Hz light source is used.
Default	Restores the default settings.

MEASUREMENT TOOLBAR

To open the Measurement Toolbar, move the mouse cursor to the top of the window. The toolbar will automatically pop up.

ICON	FUNCTION	ICON	FUNCTION
	Float/Fix switch of the Measurement Toolbar		Parallel lines
<input checked="" type="checkbox"/> Visible	Show/Hide Measurement Objects		Rectangle
Nanometer (nm) ▾	Select the desired Measurement Unit		Ellipse
4X ▾	Select Magnification for Measurement after Calibration		5 Point Ellipse: drawing up of an ellipse using five points
	Object Select		Circle: center and radius
	Angle		3 Point Circle: drawing up of a circle using three points
	4 Point Angle		Annulus
	Point		Two Circles and its Center Distance: drawing up of two circles and the measurement of the distance between the centers of these two circles
	Arbitrary Line		Two 3 Point Circles and the Center Distance: drawing up of two three-point-circles and the measurement of the distance between the centers of these circles
	3 Point Line		Arc
	Horizontal Line		Text
	Vertical Line		Polygon
	3 Point Vertical Line		Curve
	Scale Bar		Measurement Setup
	Arrow		Deletes all of the measurement objects

Executes Calibration with the help of a micrometer to determine the corresponding relation between magnification and resolution as well as between the measurement unit and the sensor pixel size.

Exits from Measurement mode

Exports the Measurement information to CSV file(*.csv)

Use these icons to select Move Left, Move Right, Move Up, Move Down, Color Adjustment, or Delete.

Select a figure to draw up; click at the desired initial point on the image, click at the desired end point. The selected figure will be drawn up on the image and the figure parameters will be displayed.


CAMERA CONTROL TOOLBAR

To open the Camera Control Toolbar, move the mouse cursor to the bottom of the window. The control toolbar will automatically pop up.

ICON	FUNCTION	ICON	FUNCTION
	Zoom In the Video Window		Zoom Out the Video Window
	Horizontal Flip		Vertical Flip
	Color/Gray		Video Freeze
	Display Cross Line		Browses the images and videos on the SD Card
	Settings		Check the Version of MAGUSView

SETTINGS

MEASUREMENT

Global	Precision	Sets the required number of digits after the decimal point. The measurement result value will be displayed with the required precision.
Calibration	Line Width	Defines width of the lines for calibration.
	Color	Defines color of the lines for calibration.
	EndPoint	Type: defines shape of the endpoints of lines for calibration: Null means no endpoints, Rectangle means rectangle shape of endpoints. It facilitates the object alignment.

Point, Angle, Arbitrary Line, Horizontal Line, Vertical Line, Rectangle, Circle, Ellipse, Annulus, Two Circles, Polygon, Curve.

You can Left-click the  and select Measurement to open the corresponding window of settings. You can specify the parameters of the required measurement objects.

MAGNIFICATION

Name	Sets the magnification item name when calibration, such as 4X, 10X, 100X, etc. Other information could be added too, for example, microscope model, operator name, etc.
Resolution	Sets the pixels per meter.
Clear All	Clears all of the calibrated magnification items.
Delete	Deletes the selected magnification item.

IMAGE FORMAT

Measurement Object Saving Method	Burn-in Mode: merges the measurement objects into one current image. Further edition of the measurement objects is impossible. The changes are irreversible.
	Layered Mode: saves the measurement objects at different layers of the current image into the target file. The measurement objects in the target file can be edited with the use of corresponding software on the PC. The changes are reversible.

VIDEO

Fast Forward/Reverse Interval	Sets the time interval of the playback of video files.
Video Encode	You can select the video encoding format H264 or H265.

STORAGE

Storage Device	SD Card is the only supported storage device.
File System Format of the Storage Device	You can select the file system format: FAT32, exFAT or NTFS. Use PC to format the SD cards and switch between the file system formats. The maximum video file size of a single file is 4GB. Unknown Status: SD card not detected or the file system is not identified.


FILES

Image File Name Auto: Saves the image file automatically with the specified prefix.
Manual: You need to name the image file before saving.

Video File Name Auto: Saves the video file automatically with the specified prefix.
Manual: You need to name the video file before saving.

The maximum video file size is 4GB. Multiple video files will be automatically generated in case of long video recording.

MISCELLANEOUS

Auto Exposure Sets the maximum exposure time during the auto exposure process. You can set a lower exposure time value to get a faster frame rate during auto exposure.

ROI Color Selects a line color of the rectangle Region of Interest (ROI).

Camera Parameters Import Imports the camera parameters from the SD card to use the previously exported camera parameters.

Camera Parameters Export Exports the camera parameters to the SD card to use the previously exported camera parameters.

Reset to factory defaults Restores the camera parameters to the default settings.

LANGUAGE

You can also select a language from the list.

THE KIT INCLUDES

MAGUS microscope digital camera, 12V, 1A power adapter, HDMI cable, USB mouse, 32GB SD card, mounting plate and screws, user manual, and a warranty card.

SPECIFICATIONS

MAGUS CHD10	
Maximum resolution	1920x1080px
Megapixels	2
Sensor	1/2.8" (5.57x3.13mm) Sony CMOS sensor
Pixel size	2.9x2.9µm
Light sensitivity	1300mv with 1/30s
Exposure	0.01ms–1000ms
Video recording	+
Frame rate, fps at resolution	60@1920x1080px
Image format	*.jpeg
Video format	*.mp4
Shutter type	ERS
Interface	HDMI
Body	metal
Operating temperature range	–10...50°C

The manufacturer reserves the right to make changes to the product range and specifications without prior notice.

CARE AND MAINTENANCE

- Do not disassemble the camera on your own.
- Keep the camera away from moisture; do not use it in rainy weather.
- Keep the camera away from shock or excessive pressure.
- Do not overtighten the locking screws.
- Keep the camera away from hazardous environment, home and car heaters, incandescent lamps or open fire.
- When cleaning the lenses, first blow any dust or debris off the surface or wipe them off with a soft brush. Then wipe the lens with a soft cleaning cloth slightly damp with spirits or ether.
- Seek medical advice immediately if a small part or a battery is swallowed.

MAGUS WARRANTY

Magus optic products carry a **5-year warranty** against defects in materials and workmanship. All Magus accessories are warranted to be free of defects in materials and workmanship for **six months** from the purchase date. The warranty entitles you to the free repair or replacement of the Magus product in any country where a Levenhuk office is located if all the warranty conditions are met.

For further details, please visit: www.magusmicro.com

If warranty problems arise, or if you need assistance in using your product, contact the local Levenhuk branch.

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MAGUS® is a registered trademark of Levenhuk, Inc.
© 2006–2024 Levenhuk, Inc. All rights reserved.